

Unravelling the Complexity of Fragrances in HPC* Products Using

TruTOF GC-TOF-MS with Automated Data Deconvolution

Nick Jones,¹ Marcel Kneepkens,¹ Izabela Malyszko,¹ Sjaak de Koning,² and Lorraine Kay,²

¹Global Fragrance Centre, Unilever Research and Development, Port Sunlight, UK,

²LECO European Technical Centre (ETC), Separation Science, Mönchengladbach, Germany.

Introduction

The composition and analysis of fragrance components in home and personal care (HPC) products is very complex and unquestionably time consuming.

Traditionally, significant GC-MS time may be consumed by running multiple separations using different phases. The resulting chromatograms are subject to labour intensive data interpretation, using retention index and mass spectral libraries, to accurately resolve and quantify the composition of perfume raw materials, including essential oils in the products.

Now, a new analytical approach has been evaluated. Using a LECO TruTOF GC-TOF-MS, with automated true signal data deconvolution, it is possible to significantly increase the throughput of such analyses using single separations. Additionally, in some cases where either there are a lot of interferences from the matrix or a complex mixture of essential oils, the ability to deconvolute full range mass spectra provides extra information. This will also increase our ability to accurately characterize and interpret our data to further unravel the complexity of fragrances in such

products.

A skin care lotion, containing a relatively complex perfume including a number of essential oils was chosen to assess the ability of LECO TruTOF GC-TOF-MS. The goal was to increase both the throughput and quality of output, compared to the conventional multi-run analysis approach described. The skin care lotion analysed, was particularly appropriate, because of the presence of a high level of alkane (petroleum jelly) components. Such materials are significant interferences in fragrance analysis, as a result of their similarity in polarity and volatility to a high number of middle and base note perfume components.

Sample Selection and Introduction

A solvent extraction and proprietary clean-up procedure was used to extract perfume components from the skin lotion. Internal standards were added prior to the extraction/clean-up process for quantification purposes. The resulting sample was then concentrated by Vigreux distillation for analysis. Analysis of the sample was then performed by the traditional multi-phase/run approach described, using

single quadrupole GC-MS systems, as well as by a single separation using a LECO TruTOF GC-TOF-MS system.

Results

Analysis of data collected using the existing 2-separation, single quadrupole MS measurements, processed manually, was compared to data collected from a single sample injection and separation using the TruTOF GC-TOF-MS system. This data was automatically deconvoluted using LECO's ChromaTOF software and incorporated Unilever R&D, Global Fragrance Centre mass spectral fragrance libraries.

The output of the TruTOF analysis was compared to that of the standard, more time consuming fragrance analysis, which as a fully validated and well used procedure, provides an excellent benchmark. Key outputs from the analysis of fragrances in HPC products are inclusion levels, top, middle and base notes, essential oils, musks and approximate costs (as shown in Table 1). Therefore, measures of these were chosen for the comparison. The single injection TruTOF approach compared very favourably with the standard approach. In addition to this, a number of additional musk and essential oil

Table 1:

Analysis	TruTOF	Current output
Inclusion level (%)	0.16	0.16
Notes		
Top	35.9	32.1
Middle	51.9	63.1
Bottom	8.5	4.8
Essential oils		
Orange oil	1.6	1.6
Ylang ylang	1.5	1.42
Musks		
Polycyclic	6.4	4.37
Macrocyclic	0.2	ND
Cost (Euro/kg)	13.3	12.5

Figure 1: TIC showing the 'hump' of interfering alkane components, eluting from 4000 to 6000 seconds, which without data deconvolution of full range mass spectra, mask a number of fragrance and essential oil components, rendering data analysis very time consuming and challenging.





Figure 2: Deconvoluted chromatograms revealing a number of peaks eluting underneath and closely alongside the principle peak (333) of Germacrene D, a component of ylang ylang essential oil. Peak 332 is a component of the synthetic, Methyl Ionone Alpha.



components were identified underneath the 'hump' of alkane interferences.

The TIC trace in Figure 1 does not lend itself to easy interpretation. However, Figure 2 shows a zoom of the mass traces that are hidden under the TIC and seen following the automated deconvolution data processing methods. The deconvoluted chromatograms also reveal a number of peaks eluting underneath and closely alongside the principle peak (333) of Germacrene D, a component of ylang ylang essential oil. Peak 332 is a component of the synthetic, methyl ionone alpha and is represented by the pink mass trace below.

Conclusion

This study suggests that the use of a LECO TruTOF GC-TOF-MS system coupled with ChromaTOF deconvolution software, may provide an accurate and higher throughput alternative to other, multi-phase, multi-run separations that require significantly more resources, both in terms of man hours for data analysis and instrument operation time.

In addition, the ability to maintain full sensitivity over the full mass range and the ability to collect spectra at up to 80 spectra/second, yielding effective data deconvolution, may enable a richer and more informative data set and, to a great extent, eliminate the risk of losing important information during the analysis of HPC products that contain a high level of interferences, such as the skin lotion analysed here.



Delivering the Right Results

LECO European Technical Centre (ETC)

Marie-Bernays-Ring 31, 41199
Mönchengladbach, Germany
tel. +49 2166 687 104 or 107
fax +49 2166 687 148
E-mail: SepSci@leco.de
Website: www.leco.com

Article Reprinted from ©Dec 17, 2009 issue of

The **Column**

the interactive e-publication for the european separation science industry
AN ADVANSTAR PUBLICATION