

Non-Targeted Analysis of Complex Environmental Samples Using Thermal Desorption, Multidimensional Gas Chromatography, Multi-Mode Ionization Methods, and High-Resolution Time-of-Flight Mass Spectrometry

David E. Alonso, John Hayes, Liz Humston-Fulmer, and Joseph E. Binkley | LECO Corporation, Saint Joseph, Michigan, USA




Introduction

- Confident analysis of environmental samples (e.g., air, water, dust/soil) is critical for the monitoring of harmful substances such as persistent organic pollutants
- Many of these harmful chemical substances have been classified as carcinogens and endocrine disruptors
- Study Objective: To use GCxGC-HRTOFMS for the analysis of complex environmental samples.

Challenge: Complex environmental samples are difficult to analyze

- Very large number of sample constituents
- Wide concentration range
- Varying physiochemical properties

Solution: The utilization of powerful analytical discovery tools



Instrument Specifications




Figure 1. Pegasus® HRT® 4D and Multimode Source.

Table 1: Instrument Parameters

Gas Chromatograph	Agilent 7890B with LECO Dual Stage QuadJet™ Modulator
Sample Introduction A	TD, Split 10:1 (3 mins) → Splitless, 50 °C to 300 °C at 10 °C/s
Sample Introduction B	Pyr, 10:1 → Splitless (3 mins), 50 °C to 600 °C at 60 °C/s
Carrier Gas	He @ 1.4 mL/min, Corrected Constant Flow
Primary Column	Rxi-5ms, 30 m x 0.25 mm i.d. x 0.25 µm
Secondary Column	Rxi-17sil ms, 1.3 m x 0.10 mm x 0.10 µm
Temperature Program	50 °C (2 min) ramp 10 °C/min to 250 °C, then ramp to 300 °C at 20 °C/min to 300 °C (Hold 15 mins) Secondary oven maintained +5 °C relative to primary oven
Modulation Period	3.0 seconds; modulator maintained +15 °C relative to secondary oven
Transfer Line	300 °C
Mass Spectrometer	LECO Pegasus HRT® 4D
Source Temperature	EI, 250 °C; PCI, ECNI, 165 °C
Acquisition Mode	High Resolution; Re ≥ 25,000 for m/z 219, Mass Accuracy ≤ 1 ppm
Ionization	EI, PCI/ECNI (Reagent Gas = CH ₄)
Mass Range (m/z)	EI 50-1000; PCI 60-1000; ECNI 30-1000
Acquisition Rate	200 sps

Monitoring where we live and work

1) Dust




Figure 2. Contour plot of a complex dust sample.




Figure 3. Spectral Analysis Tools (SAT) for complex sample simplification. Mass defect plot of dust sample.




Figure 4. Filtered (XICP) contour plot displaying a homologous series of alkylmethylbenzylamines.

Monitoring the food we consume

2) Food (Fish)




Figure 5. Contour plot displaying polychlorinated biphenyl region and selected halogenated organic compounds.

Table 2: Selected Halogenated Compounds

Name	R.T. (s)	Formula	PPM	Similarity	Name	R.T. (s)	Formula	PPM	Similarity
Benzene, hexachloro	1595.2860	C ₆ Cl ₆	0.87	928	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2340.2400	C ₁₂ H ₂ Cl ₈	0.12	717
Phenol, 2-chloro-4-methoxy	1000.3800	C ₆ H ₅ ClO ₂	1.37	932	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2340.2400	C ₁₂ H ₂ Cl ₈	0.12	807
2,2',3'-Trichloro-1,1'-biphenyl	1915.2385	C ₁₂ H ₃ Cl ₃	0.32	787	2,3,3',4,4'-Pentachloro-1,1'-biphenyl	2345.2462	C ₁₂ H ₂ Cl ₆	0.33	872
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	1970.2399	C ₁₂ H ₂ Cl ₈	1.55	795	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2350.2462	C ₁₂ H ₂ Cl ₈	0.33	872
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	1995.2352	C ₁₂ H ₂ Cl ₈	1.06	679	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2355.2462	C ₁₂ H ₂ Cl ₈	0.33	872
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2030.2420	C ₁₂ H ₂ Cl ₈	0.76	686	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2360.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2045.2440	C ₁₂ H ₂ Cl ₈	0.39	942	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2365.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2065.2420	C ₁₂ H ₂ Cl ₈	0.57	933	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2370.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2105.2420	C ₁₂ H ₂ Cl ₈	0.99	726	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2375.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2150.2420	C ₁₂ H ₂ Cl ₈	0.08	905	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2380.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2155.2420	C ₁₂ H ₂ Cl ₈	0.68	879	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2385.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2160.2420	C ₁₂ H ₂ Cl ₈	0.75	931	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2390.2462	C ₁₂ H ₂ Cl ₈	0.33	931
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2175.2420	C ₁₂ H ₂ Cl ₈	0.08	905	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2410.2462	C ₁₂ H ₂ Cl ₈	0.33	884
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2175.2420	C ₁₂ H ₂ Cl ₈	0.08	905	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2415.2462	C ₁₂ H ₂ Cl ₈	0.33	884
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2175.2420	C ₁₂ H ₂ Cl ₈	0.61	911	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2420.2462	C ₁₂ H ₂ Cl ₈	0.33	884
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2185.2395	C ₁₂ H ₂ Cl ₈	0.81	922	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2475.2320	C ₁₂ H ₂ Cl ₈	0.03	917
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2215.2375	C ₁₂ H ₂ Cl ₈	0.18	901	2,2',3,3',4,4',5,5'-Hexachloro-1,1'-biphenyl	2495.2345	C ₁₂ H ₂ Cl ₆	0.79	908
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2260.2320	C ₁₂ H ₂ Cl ₈	1.17	901	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2495.2345	C ₁₂ H ₂ Cl ₈	0.39	890
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2265.2340	C ₁₂ H ₂ Cl ₈	1.32	780	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2505.2340	C ₁₂ H ₂ Cl ₈	0.14	849
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2275.2340	C ₁₂ H ₂ Cl ₈	0.87	873	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2515.2340	C ₁₂ H ₂ Cl ₈	0.27	926
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2275.2340	C ₁₂ H ₂ Cl ₈	0.89	889	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2525.2340	C ₁₂ H ₂ Cl ₈	0.83	934
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2290.2385	C ₁₂ H ₂ Cl ₈	0.43	900	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2535.2340	C ₁₂ H ₂ Cl ₈	0.53	924
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2305.2358	C ₁₂ H ₂ Cl ₈	N/A	780	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2545.2355	C ₁₂ H ₂ Cl ₈	0.46	926
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2310.2360	C ₁₂ H ₂ Cl ₈	0.5	942	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2560.2399	C ₁₂ H ₂ Cl ₈	0.18	885
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2320.2385	C ₁₂ H ₂ Cl ₈	0.32	820	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2575.2360	C ₁₂ H ₂ Cl ₈	0.12	925
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2335.2355	C ₁₂ H ₂ Cl ₈	0.58	938	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2605.2340	C ₁₂ H ₂ Cl ₈	0.22	925
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2340.2340	C ₁₂ H ₂ Cl ₈	0.12	703	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2610.2340	C ₁₂ H ₂ Cl ₈	0.31	852
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	2340.2395	C ₁₂ H ₂ Cl ₈	0.35	903					

Complex Environmental Samples

Monitoring the air we breathe

3) Air Particulate Matter



Air Samples: Collection and Introduction




Figure 6. Active air sampler (left). Placed in urban and rural sites (24-hour sampling). GL Sciences Optic-4 Inlet (right) used for Thermal Desorption (TD) and Pyrolysis (Pyr).

Air Sample A (TD)




Figure 7. Contour plot displaying areas of group clustering for an air sample.

Table 3: Representative Compounds in Sample A

Name	Formula	R.T. (s)	Similarity	PPM	Name	Formula	R.T. (s)	Similarity	PPM	Name	Formula	R.T. (s)	Similarity	PPM
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	1.77	Benzoic acid, 2-methyl-	C ₈ H ₇ O ₂	573.1445	843	0.8	2-Ethyl-1-dodecanol	C ₁₄ H ₂₈ O	999.2271	833	N/A
Chloro-2,3-dibiphenyl	C ₁₂ H ₉ Cl	195.1210	824	-0.26	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	195.1210	824	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	824	0.25
Pyridine	C ₅ H ₅ N	153.1408	830	0.77	2-Propenoic acid, 1-methyl-, (E)	C ₅ H ₇ O ₂	588.1336	847	1.29	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	802	N/A
Pyridine	C ₅ H ₅ N	153.1408	793	0.25	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	793	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	N/A
Cyclopentane, 1,2,4-trimethyl-	C ₅ H ₈	241.2016	791	0.99	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	843	1.38	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	1.42
1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	871	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	942	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	1.42
Maleic anhydride	C ₄ H ₂ O ₃	309.1064	938	-0.24	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	0.24	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	802	-0.4
4-Cyclopentene-1,2-dione	C ₅ H ₄ O ₂	333.1228	892	-0.88	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	873	-0.54	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	848	0.81
2-Pyridone	C ₅ H ₅ NO	128.1286	793	0.25	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	827	0.46	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	841	-0.15
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	831	0.74	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethenyl-	C ₅ H ₅ N	181.2740	893	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	837	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Benzenethiol	C ₆ H ₆ S	106.1286	932	0.36	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	836	0.77	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	841	0.82
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	901	0.68	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	937	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.49
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	0.76	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.03
1-Hexanol, 2-ethyl-	C ₈ H ₁₈ O	165.1408	813	0.87	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	930	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
2-Pyridone	C ₅ H ₅ NO	128.1286	813	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	931	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	901	0.68	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	937	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.49
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	0.76	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.03
1-Hexanol, 2-ethyl-	C ₈ H ₁₈ O	165.1408	813	0.87	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	930	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
2-Pyridone	C ₅ H ₅ NO	128.1286	813	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	931	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	901	0.68	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	937	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.49
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	0.76	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.03
1-Hexanol, 2-ethyl-	C ₈ H ₁₈ O	165.1408	813	0.87	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	930	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
2-Pyridone	C ₅ H ₅ NO	128.1286	813	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	931	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	901	0.68	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	937	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.49
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	0.76	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.03
1-Hexanol, 2-ethyl-	C ₈ H ₁₈ O	165.1408	813	0.87	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	930	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
2-Pyridone	C ₅ H ₅ NO	128.1286	813	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	931	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	901	0.68	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	937	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.49
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	0.76	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.03
1-Hexanol, 2-ethyl-	C ₈ H ₁₈ O	165.1408	813	0.87	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	930	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
2-Pyridone	C ₅ H ₅ NO	128.1286	813	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	931	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	901	0.68	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	937	0.4	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.49
Hexachlorocyclopentadiene	C ₅ Cl ₈	141.2444	893	0.76	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.03
1-Hexanol, 2-ethyl-	C ₈ H ₁₈ O	165.1408	813	0.87	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	930	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
2-Pyridone	C ₅ H ₅ NO	128.1286	813	0.71	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	808	1.01
1-Cyclopentene-1,1-dione	C ₅ H ₄ O ₂	339.1061	891	0.51	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	931	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	807	0.58
Pyridine, 2-ethyl-	C ₅ H ₅ N	181.2740	905	0.84	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	938	N/A	1,1'-Biphenyl, 2,2',3,3',4,4',5,5'-octachloro	C ₁₂ H ₂ Cl ₈	1995.2352	809	0.94