



Delivering the Right Results

For Immediate Release
October 10, 2007

Time-of-Flight Capabilities For Any Laboratory and Budget

St. Joseph, MI—The TruTOF™ HT TOFMS offers all the power, speed, and software features of our Pegasus® HT in a cost-competitive benchtop design. Continuous full-range mass acquisition rates up to 80 spectra/second give you accurate identification of unknown and targeted compounds in complex samples. Exclusive ChromaTOF® software with True Signal Deconvolution™ provides fast and easy data analysis and clean library searches. The TruTOF also includes environmental reporting software, and both electron and chemical ionization sources for increased laboratory productivity.

For more information contact:

LECO Corporation • 3000 Lakeview Avenue • St. Joseph, MI 49085 • U.S.A.
Phone: (800) 292-6141 • Fax: (269) 982-8977 • email: info@leco.com • www.leco.com

LECO is a registered trademark of LECO Corporation.