

# MSX250 A-Series

## Specification Sheet

These versatile benchtop sectioning machines offer both linear and radial feed configurations in efficient benchtop models.

Supports up to 12 in. (300 mm) blade @ 5.4 HP (4 kW) with integrated electronic braking system

Sealed spindle design for optimal speed, precise sections, and trouble-free operation

Easily accessible sealed electronic enclosure



Radial arm with manual locking feature for table sectioning



Channeled coolant return for high flow rates and minimal debris build-up




Convenient easy-to-operate control panel

Capable of storing up to 10 applications

Automatic pulse feed rate with load control maximizes sectioning efficiency without burning samples

### MSX250A1



13.8 in. (350 mm) x 9.8 in. (250 mm) T-slotted table

**LECO**<sup>®</sup>  
Delivering the Right Results