

MSX255M-Series

Specification Sheet

These versatile benchtop sectioning machines offer both linear and radial feed configurations in efficient benchtop models.

Sealed spindle design for optimal speed, precise sections, and trouble-free operation

Safety shut-off switch terminates operation when hood is opened during use


Easily accessible sealed electronic enclosure

Supports up to 12 inch (300 mm) blade @ 5.4 HP (4.0 kW) with integrated electronic braking system

Optional Table Vises



Convenient easy-to-operate control panel



T-slotted table
13.8 in. (350 mm) x
9.8 in. (250 mm)



Channeled coolant return for high flow rates and minimal debris build-up

MSX255M

LECO[®]
Delivering the Right Results