



# Selective Detection of Halogenated Pesticides in Citrus Oils with GCxGC-TOFMS and Scripts

Life Science and Chemical Analysis Centre  
LECO Corporation, Saint Joseph, Michigan USA



**Columns:**  
50 m x 0.2 mm x 0.5  $\mu$ m DB-PIONA  
2 m x 0.1 mm x 0.1  $\mu$ m BPx-50

**TOFMS:**  
35 to 350 m/z @ 100 spectra/s

For questions on this analysis  
e-mail us at: [life\\_science@leco.com](mailto:life_science@leco.com)

For a complete listing of snapshots  
and application notes visit us on the  
web at [www.leco.com](http://www.leco.com)

