



Analysis of Pump Diesel by GCxGC-TOFMS

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



Columns:

27 m x 0.25 mm x 0.25 μm Rtx-5
1 m x 0.10 mm x 0.20 μm Rtx-Wax

TOFMS:

35 to 400 m/z @ 100 spectra/s

For questions on this analysis
e-mail us at: life_science@leco.com

For a complete listing of snapshots
and application notes visit us on the
web at www.leco.com



Application Snapshot

