


Analysis of Tobacco Oil by GC-TOFMS

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



Total Ion Chromatogram (TIC) for the analysis of a tobacco oil sample. The insets show detail in crowded regions of the chromatogram.

Column:

9.4 m x 0.18 mm x 0.2 μ m DB-5

TOFMS:

35 to 300 m/z @ 20 spectra/s

For questions on this analysis e-mail us at: life_science@leco.com

For a complete listing of snapshots and application notes visit us on the web at www.leco.com

