


Note: Part numbers and standards' values may change. Consult LECO for the latest information.



GENERAL COMBUSTION SUPPLIES

766-036 Air Regulator
 502-091 Calcium Oxalate Standard (50 g)
 621-331 Ceramic Crucible—Large (5 pkg)
 529-048 Ceramic Crucible Cover—Large (5 pkg)
 502-335 Crucible Foil Liner (500/pkg)
 761-929 Crucible Tongs
 614-442-110 Crucible Turntable—Large (Kit)
 776-978 Glass Scoop
 764-216 Helium Regulator
 501-291 Oxygen Regulator
 778-898 Crucible Cover Forceps

REFERENCE MATERIALS**COAL PROX-PLUS (50 grams)**

Part No.	Desc.	% Volatile	% Ash
502-442	Low	36.0	7.2
502-443	Medium	43.5	9.7
502-444	High	37.2	7.4

COKE PROX-PLUS (50 grams)

Part No.	Desc.	% Volatile	% Ash
502-452	Prox Plus	11.7	0.37

CALCIUM OXALATE (50 grams)

Part No.	% Weight Loss		
	@ 200°C	@ 400°C	@ 850°C
502-091	12	19	30

*NOTE: Part numbers and standards' values may change.
 Consult LECO for latest information.*


United States:

Ph: 1-800-292-6141
 Fax: 269-982-8977

Canada:

Ph: 1-800-268-0215
 Fax: 905-564-6582

represented by:

**Other Countries:**

Ph: 269-983-5531
 Fax: 269-983-3850

© 2010 LECO Corporation